



Electromagnets

Technical Specifications

- Full Scale Field : ± 1 Tesla / ± 2 Tesla / ± 2.5 Tesla / ± 3.2 Tesla
- Water cooled from coil faces
- Adjustable Gap: 0-100mm
- Various Pole Face dimensions: 25mm / 50mm / 75mm / 100mm
- Excellent Uniformity: ± 0.5 % in XYZ
- Water flow and temperature sensors in each coil
- Closed Loop Operation with Hall Probe Gaussmeter
- Air cooled Bipolar Power Supply : 3kVA / 5kVA / 10kVA
- Universal USB 2.0 Interface

